



**COMPLETE PROTECTION KIT  
AGAINST COVID-19  
& CONTAMINATED AIRBORNE PARTICLES**

# 3 PLY FACE MASKS

3 Ply Face Mask helps reduce the spread of viruses, germs and bacteria. It provides relief and allows freer breathing for allergy sufferers. The 3 Ply face mask contains a high level static filtration non-woven middle layer to effectively block micro dust and bacteria.

Quality crafted for excellent protection and a comfortable fit, this face mask is ideal for a wide range of applications and wherever fluid exposure is probable. Specialized nonwoven material structure features a three-ply, pleated construction. This triple layer medical mask is a disposable mask, fluid-resistant, providing protection from droplets of infectious material emitted during coughing/sneezing/talking.



# FFP2 RESPIRATORY MASK

This Filtering half mask protects against particles FFP2 NR, without exhalation valve, limited to single shift use only. The filtering half mask is foldable, designed with a nose slide in aluminium and two adjustable head harnesses in elastic strap. The filtering media is made of melt blown reference RMF20266/2, covered by a spun bond layer reference RMF-40266 and an inner layer reference RMFLP160250.



Disposable Dust Mask, for protection against dust particles and fumes

Tested and certified to EN standard for FFP2 level

Weld fringe seal to avoid open layer around the edges

The material used offers high filtering efficiency.

Contour design ensures the compatibility of glasses/goggles and reduces fogging

# 3M™ PARTICULATE RESPIRATOR 8210, N95 MASK

Reliable, value-priced protection against a wide range of non-oily particles and engineered to provide 95% filtration efficiency. Advanced Electrostatic Media technology enhances particulate capture while reducing breathing resistance and exposure to certain airborne particles. This disposable respirator is designed for protection against particles arising from grinding, sanding, sweeping, sawing, bagging or other dusty operations. It also helps reduce inhalation exposures to certain airborne biological particles (examples: mold, Bacillus anthracis, Mycobacterium tuberculosis).



Two-strap design with welded dual point attachment provides greater durability and also helps provide a secure seal

Cushioning soft nose foam for more comfort and adjustable noseclip

Being lightweight, promotes greater worker acceptance and wear time

Compatible with a variety of protective eyewear and hearing protection

**Note:** Product images shown for representative purposes. Actual images with specifications would be provided as per the product availability at the time of placing order. We can provide re-packaged kit with selected items with your branding as part of your CSR obligation.

# N95 RESPIRATOR MASKS V-4400 SERIES



V-4400 series of respirators has superior micro fine filter media technology which protects against respirable suspended particulate matter. These respirators consist of multiple layers viz pre-filter, fine filter and skin comfort layer and should be used to protect the wearer from solid and non-volatile liquid particles and against mechanically generated particulates such as those generated from sanding, grinding, drilling, sweeping etc. These respirators have filters that can help reduce breathing in pathogenic biological airborne particulates such as influenza virus.



Filters > 95% non-oil based particles & at least > 80% solid air borne particles

Embedded nose clip – no metal exposure or loose parts

Unique fit adjusters provides optimum protection and comfort

Superior technology micro fine filter media

Latex free knitted textile elastic - longer life, skin friendly and does not deform in high temperature

## HAND SANITIZERS

Health-specialists have been emphasising upon the importance of keeping your hands germ-free, either through the use of soaps or hand sanitisers. Although meticulous hand washing is your best bet to get rid of the grime that clings to the hand, that might harbour nasty pathogens, including the coronavirus, however in situations where washing your hands is not possible, hand sanitizer can prove to be useful in eliminating dangerous microorganisms.

At present, alcohol-based sanitizers are the only known means for rapidly and effectively inactivating a wide array of potentially harmful microorganisms on hands.

The hands must be thoroughly scrubbed, making sure that the back, the place between the fingers and the underside of nails are not left out.



Available  
in 100ml,  
200ml, 500ml,  
1 litre & 5 litres  
pack

BUY AT GOVT  
PRESCRIBED PRICES

## PERSONAL PROTECTION EQUIPMENT (PPE) KIT

These are specific outfits or tools worn out by health care workers to get protected by various infectious elements. These premium quality Personal protective equipment works as a barrier between healthcare specialist and suspicious infected personal.

PPEs are not alternative to basic preventive public health measures such as hand hygiene, respiratory etiquettes. Always (if possible) maintain a distance of at least 1 meter from contacts/suspect/confirmed COVID-19 cases and always follow the laid down protocol for disposing off PPEs as detailed in infection prevention and control guideline available on website of MoHFW.



### Contents of PPE Kit

Disposable Coverall | Goggles | 3 Ply Face Mask | Nitrile Gloves  
Shoe Cover | Bio Hazard Bag

**Note:** Product images shown for representative purposes. Actual images with specifications would be provided as per the product availability at the time of placing order. We can provide re-packaged kit with selected items with your branding as part of your CSR obligation.



# COVID 19 ANTIBODY TEST KIT

The antibody tests, also known as serological tests, look for the presence of a virus through a blood test rather than a nasal or throat swab. This is usually done as an academic exercise to check whether an individual had the presence of the disease and recovered on their own, helping researchers map the spread of the virus as well as the number of people in the population who have immunity against the virus.

It may be noted that this is not a confirmatory test. The serological test is for the purpose of surveillance to generate data and understand whether people got exposed to the virus.




# INFRARED FOREHEAD THERMOMETER

This digital thermometer uses infrared energy radiated by the body to measure your temperature when 3-5 cm away from your forehead. The infrared thermometers are used to detect the body temperature in case of fever, one of the symptoms of coronavirus. Non-contact measuring is helpful to avoid cross-infection.

SAFE FAST CONVENIENT

## Specifications

Measuring Distance	0-5 cm
Display	LCD
Resolution	0.1°C/0.1°F
Measuring Range	0°-100°C
Accuracy	±0.3°C
Operating Temperature & Humidity	5°C - 40°C, 85%RH
Operating Temperature & Humidity	-20°C-55°C, 95%RH
Power Supply	2x1.5V AAA Batteries
Service Life	5 Years



# FEVER SCREENING METAL DETECTOR DOOR SOLUTION

This is a Metal detector door and thermal camera all-in-one solution. It gives real-time display of statistics extracting the temperature of forehead and detecting faces by AI algorithm through thermal imaging from 30°C to 45°C with ±0.5 °C accuracy.

## Advantages

- Hygienic touch-free fast fever screening
- Reducing human resources cost
- Preventing people with abnormal temperature entering
- Long distance verification & measurement
- Convenient information display on LCD
- Suitable for various indoor environments – Malls, Cinemas, Hospitals, Airports, Educational Institutions



**Note:** Product images shown for representative purposes. Actual images with specifications would be provided as per the product availability at the time of placing order. We can provide re-packaged kit with selected items with your branding as part of your CSR obligation.

# VENTILATORS

According to the WHO, One in six COVID-19 patients can develop breathing difficulties, potentially needing access to a ventilator. Ventilators - mechanical breathing devices - are crucial in the fight to save patients whose lungs are assailed by the virus.

## AK-8800 ICU Ventilator



An ICU Ventilator with Gas-Electricity separation ensuring safe operation with the facility of dis-assembling & disinfecting Exhalation valve to avoid cross-infection.

## AV-2000B, Ventilator



A multi-functional Ventilator applicable for children and adult in sickroom, ICU or emergency.

**Honeywell**

**NORTH** Primair® 700 Series

## POWERED AIR PURIFYING RESPIRATOR (PAPR)

Workers who perform their jobs in high-risk environments face hazards on a daily basis. A powered, air-purifying respirator (PAPR) is a type of personal protective equipment used to safeguard workers against contaminated air. PAPRs consist of a respirator in the form of a hood, or full-face mask, which takes ambient air that is contaminated with one or more type of pollutant or pathogen, actively removes (filters) a sufficient proportion of these hazards, and then delivers the clean air to the user's face and/or mouth.

The PA700 protects workers from hazards so that they breathe easy, work hard, focus on the job and feel safer.



# ABOUT AHAAN HEALTH CARE

**Ahaan Healthcare (AHC)** is an entity of Acme Media Group, a 10 year 120 Member strong national Media Marketing Group headquartered in Delhi with diverse interests in Cinema Advertising, Non-Traditional Media, Out-of-home Advertising, Digital Marketing, Exhibition, Events & Conferences, Brand Activation, Travel & Tourism and Furniture & Interior. AMCPL works across various sectors – Government, Trade Associations, Industries, hospitality, Health & Pharma, FMCG, Telecom, Automobile, E-Commerce, Tourism, Fashion and Retail.

In view of the enormous health care product needs, AMCPL has diversified into the field of health care products leveraging the power of its marketing and networking skills. AHC in addition to its manufacturing base at Mandi, Himachal Pradesh, also acts as a re-selling partner to major Healthcare product suppliers by making critical genuine certified health care products available to the end-user finding it difficult to source multiple products from one roof in a market flooded with scrupulous products.



**Corporate Office:**

Acme House, 562 A, 2nd Floor, M G Road, Lado Sarai, New Delhi 110030  
**Delhi:** 95605 48282, **Chandigarh:** 8447127316, **Kolkata:** 98366 23468,  
**Bihar & Jharkhand:** 99349 83479, **North East:** 91277 53678,  
**Chennai & Bangalore:** 97394 31436.  
**E:** info@ahaanhealthcare.com, **W:** ahaanhealthcare.com